

# XE-RET

## HEMATOLOGY CONTROL

Reference Values provided by DYMIND

**CONTROL**

LOT 34091

2024-08-08

2024-11-05

Applicable Instruments	Parameter	Unit	LOT 34091001	LOT 34091002	LOT 34091003
DYMIND	RET#	$\times 10^{12}/L$	0.0352 $\pm 0.0176$	0.2237 $\pm 0.1316$	0.2883 $\pm 0.1200$
DH-600 DH-602 DH-605 DH-610 DH-612 DH-615 (Technical File Version A2.0 or higher)	RET%	%	0.64 $\pm 0.32$	4.91 $\pm 1.96$	12.27 $\pm 4.50$
	LFR	%	90.0 $\pm 27.0$	90.1 $\pm 27.0$	89.7 $\pm 27.0$
	MFR	%	7.7 $\pm 11.0$	8.1 $\pm 12.0$	8.2 $\pm 12.0$
	HFR	%	2.3 $\pm 3.0$	1.8 $\pm 3.0$	2.1 $\pm 3.0$
	IRF	%	10.0 $\pm 15.0$	9.9 $\pm 15.0$	10.3 $\pm 16.0$
	RHE	pg	19.0 $\pm 8.0$	17.0 $\pm 7.0$	16.5 $\pm 7.0$

**【NOTE】**

1. The controls should be stored in refrigerator (2°C~8°C). After opening, it will keep stable for 14 days when it is stored airtight at 2°C~8°C.
2. Before mixing and running the control after take it out from the refrigerator, please keep it at least 15 minutes until reaching room temperature(15°C~30°C).
3. Controls must be well mixed before using. Please mix gently, to avoid cells rupture and/or generating bubbles.
4. After using, put the controls back into the refrigerator to prevent contamination and evaporation.



34091001



34091002



34091003

# XE-RET

## HEMATOLOGY CONTROL

Reference Values provided by DYMIND

**CONTROL**

LOT 34091

2024-08-08

2024-11-05

Applicable Instruments	Parameter	Unit	LOT 34091001	LOT 34091002	LOT 34091003
DYMIND DM60 Vet DM61 Vet DM62 Vet DM63 Vet (Technical File Version A1.4 or higher)	RET#	$\times 10^{12}/L$	0.0113 $\pm 0.0090$	0.1753 $\pm 0.0877$	0.1527 $\pm 0.0763$
	RET%	%	0.21 $\pm 0.17$	3.94 $\pm 1.97$	6.51 $\pm 3.26$

**【NOTE】**

1. The controls should be stored in refrigerator (2°C~8°C). After opening, it will keep stable for 14 days when it is stored airtight at 2°C~8°C.
2. Before mixing and running the control after take it out from the refrigerator, please keep it at least 15 minutes until reaching room temperature(15°C~30°C).
3. Controls must be well mixed before using. Please mix gently, to avoid cells rupture and/or generating bubbles.
4. After using, put the controls back into the refrigerator to prevent contamination and evaporation.



34091001



34091002



34091003