

Instrument	Parameter	Low	Normal	High
EDAN DS-500 DS-500i DS-580 DS-580i DS-500C	WBC × 10 ⁹ /L	3.71 ± 0.50	8.57 ± 1.00	19.75 ± 2.50
	NEU# × 10 ⁹ /L	1.86 ± 0.40	4.90 ± 0.66	12.76 ± 2.00
	NEU% %	50.2 ± 9.0	57.2 ± 8.0	64.7 ± 9.0
	LYM# × 10 ⁹ /L	1.33 ± 0.40	2.36 ± 0.66	3.56 ± 1.50
	LYM %	35.8 ± 9.0	27.5 ± 8.0	18.0 ± 9.0
	MON# × 10 ⁹ /L	0.26 ± 0.14	0.54 ± 0.31	1.19 ± 1.09
	MON %	7.1 ± 4.0	6.3 ± 4.0	6.0 ± 6.0
	EOS# × 10 ⁹ /L	0.26 ± 0.16	0.77 ± 0.56	2.24 ± 1.25
	EOS %	6.9 ± 5.0	9.0 ± 5.0	11.3 ± 8.0
	BAS# × 10 ⁹ /L	1.0 ± 1.0	0.00 ± 1.0	0.00 ± 1.0
	BAS %	1.0 ± 1.0	0.0 ± 1.0	0.0 ± 1.0
	RBC × 10 ¹² /L	2.37 ± 0.24	4.66 ± 0.24	5.25 ± 0.5
	HGB g/L	58 ± 6	131 ± 6	164 ± 10
	HCT %	18.6 ± 3	41.6 ± 3	50.8 ± 5
	MCV fL	78.6 ± 5	89.1 ± 5	96.9 ± 8
	MCH pg	24.3 ± 3	28.1 ± 3	31.3 ± 3.5
	MCHC g/L	309 ± 30	315 ± 30	323 ± 30
	RDW-SD fL	51.7 ± 8	53.9 ± 8	55.6 ± 12
	RDW-CV %	16.9 ± 5	15.6 ± 5	14.6 ± 5
	PLT × 10 ⁹ /L	63 ± 20	291 ± 40	584 ± 60
PDW fL	11.0 ± 3	10.4 ± 3	10.2 ± 3	
MPV fL	9.0 ± 3	8.8 ± 3	8.7 ± 3	
PCT %	0.06 ± 0.06	0.26 ± 0.1	0.51 ± 0.3	
PLCR %	24.03 ± 15	20.93 ± 15	20.26 ± 15	
PLCC × 10 ⁹ /L	20 ± 20	61 ± 30	118 ± 40	

【NOTE】

- 1) The controls should be stored in the refrigerator (2~8°C). After opening, it will keep stable for 14 days when it is stored airtight at 2~8°C.
- 2) Please equilibrate the controls to room temperature (15~30°C) before using it.
- 3) Controls must be well mixed before use. Please mix gently to avoid cells rupturing or generating bubbles.
- 4) After using, put the controls back into the refrigerator (2~8°C) to prevent contamination and evaporation.

ASSAY VALUES AND EXPECTED RANGES

LOT EH0924

2024-11-05

Instrument	Parameter	Low	Normal	High
EDAN	WBC $\times 10^9/L$	3.47 \pm 0.60	8.16 \pm 1.00	18.72 \pm 2.50
	NEU# $\times 10^9/L$	1.55 \pm 1.0	4.33 \pm 3.80	11.75 \pm 7.60
	NEU%	46.2 \pm 10.0	53.3 \pm 9.0	62.4 \pm 8.0
H60	LYM# $\times 10^9/L$	1.20 \pm 1.0	2.09 \pm 1.5	3.18 \pm 2.0
H66	LYM%	34.7 \pm 10.0	25.6 \pm 9.0	17.0 \pm 8.0
H68	MON# $\times 10^9/L$	0.32 \pm 0.20	0.79 \pm 0.50	1.46 \pm 1.10
H69	MON%	9.3 \pm 4.0	9.7 \pm 5.0	7.8 \pm 6.0
H60S	EOS# $\times 10^9/L$	0.30 \pm 0.20	0.85 \pm 0.50	2.23 \pm 1.30
H66S	EOS%	8.8 \pm 5.0	10.4 \pm 6.0	11.9 \pm 7.0
H68S	BAS# $\times 10^9/L$	0.10 \pm 0.10	0.10 \pm 0.10	0.10 \pm 0.10
H69S	BAS%	1.0 \pm 1.0	1.0 \pm 1.0	1.0 \pm 1.0
	RBC $\times 10^{12}/L$	2.28 \pm 0.30	4.51 \pm 0.40	5.10 \pm 0.50
	HGB g/L	59 \pm 6	131 \pm 8	164 \pm 8
	HCT%	18.0 \pm 3.0	40.3 \pm 4.0	49.5 \pm 5.0
	MCV fL	78.9 \pm 5.0	89.4 \pm 5.0	97.1 \pm 6.0
	MCH pg	25.9 \pm 2.5	29.0 \pm 2.5	32.1 \pm 2.5
	MCHC g/L	329 \pm 30	325 \pm 30	331 \pm 30
	RDW-CV%	16.8 \pm 3.0	15.7 \pm 3.0	14.8 \pm 3.0
	RDW-SD fL	51 \pm 10	53 \pm 10	54 \pm 12
	PLT $\times 10^9/L$	56 \pm 30	279 \pm 45	507 \pm 70
	MPV fL	11.2 \pm 3.0	10.5 \pm 3.0	10.4 \pm 3.0
	PDW fL	14.2 \pm 4.0	13.0 \pm 4.0	12.9 \pm 4.0
	PCT%	0.06 \pm 0.06	0.29 \pm 0.20	0.53 \pm 0.20
	PLCR%	34.2 \pm 10.0	29.0 \pm 10.0	28.4 \pm 10.0
	PLCC $\times 10^9/L$	19 \pm 19	81 \pm 25	144 \pm 40

【NOTE】

- 1) The controls should be stored in the refrigerator (2~8°C). After opening, it will keep stable for 14 days when it is stored airtight at 2~8°C.
- 2) Please equilibrate the controls to room temperature (15~30°C) before using it.
- 3) Controls must be well mixed before use. Please mix gently to avoid cells rupturing or generating bubbles.
- 4) After using, put the controls back into the refrigerator (2~8°C) to prevent contamination and evaporation.

ASSAY VALUES AND EXPECTED RANGES

LOT EH0924

2024-11-05

Instrument	Parameter	Low	Normal	High
EDAN	WBC $\times 10^9/L$	3.61 \pm 0.60	8.38 \pm 1.00	19.28 \pm 2.50
	NEU# $\times 10^9/L$	1.62 \pm 1.0	4.10 \pm 3.80	11.69 \pm 7.60
H60 CRP	NEU%	46.3 \pm 10.0	48.9 \pm 9.0	60.1 \pm 8.0
H66 CRP	LYM# $\times 10^9/L$	1.27 \pm 1.0	2.37 \pm 1.5	3.57 \pm 2.0
H60	LYM%	35.4 \pm 10.0	28.3 \pm 9.0	18.5 \pm 8.0
CRP&SAA	MON# $\times 10^9/L$	0.35 \pm 0.30	0.85 \pm 0.50	1.71 \pm 1.30
H66	MON%	9.8 \pm 4.0	10.2 \pm 5.0	8.9 \pm 6.0
CRP&SAA	EOS# $\times 10^9/L$	0.27 \pm 0.20	0.97 \pm 0.70	2.21 \pm 1.30
H60S CRP	EOS%	7.6 \pm 5.0	11.6 \pm 6.0	11.5 \pm 7.0
H66S CRP	BAS# $\times 10^9/L$	0.10 \pm 0.10	0.10 \pm 0.10	0.10 \pm 0.10
H60S	BAS%	1.0 \pm 1.0	1.0 \pm 1.0	1.0 \pm 1.0
CRP&SAA	RBC $\times 10^{12}/L$	2.07 \pm 0.30	4.55 \pm 0.40	5.10 \pm 0.50
H66S	HGB g/L	59 \pm 6	131 \pm 8	165 \pm 8
CRP&SAA	HCT%	16.4 \pm 3.0	40.7 \pm 4.0	49.1 \pm 5.0
H80	MCV fL	79.2 \pm 5.0	89.5 \pm 5.0	96.3 \pm 6.0
H88	MCH pg	28.3 \pm 2.5	29.4 \pm 2.5	32.4 \pm 2.5
H80S	MCHC g/L	357 \pm 30	321 \pm 30	336 \pm 30
H88S	RDW-CV%	17.3 \pm 3.0	16.4 \pm 3.0	15.8 \pm 3.0
	RDW-SD fL	52 \pm 10	54 \pm 10	57 \pm 12
	PLT $\times 10^9/L$	58 \pm 30	278 \pm 45	515 \pm 70
	MPV fL	10.2 \pm 3.0	10.8 \pm 3.0	10.7 \pm 3.0
	PDW fL	10.4 \pm 4.0	13.0 \pm 4.0	12.6 \pm 4.0
	PCT%	0.06 \pm 0.06	0.30 \pm 0.20	0.55 \pm 0.20
	PLCR%	25.9 \pm 10.0	28.9 \pm 10.0	28.4 \pm 10.0
	PLCC $\times 10^9/L$	15 \pm 15	80 \pm 25	146 \pm 40

【NOTE】

- 1) The controls should be stored in the refrigerator (2~8°C). After opening, it will keep stable for 14 days when it is stored airtight at 2~8°C.
- 2) Please equilibrate the controls to room temperature (15~30°C) before using it.
- 3) Controls must be well mixed before use. Please mix gently to avoid cells rupturing or generating bubbles.
- 4) After using, put the controls back into the refrigerator (2~8°C) to prevent contamination and evaporation.

ASSAY VALUES AND EXPECTED RANGES

LOT EH0924

2024-11-05

Instrument	Parameter	Low	Normal	High
EDAN H60 Vet H66 Vet H68 Vet H69 Vet	WBC × 10 ⁹ /L	3.47 ± 0.60	8.10 ± 1.00	18.47 ± 2.50
	NEU# × 10 ⁹ /L	1.49 ± 1.0	4.21 ± 3.8	11.44 ± 7.60
	NEU% %	44.6 ± 10.0	52.2 ± 9.0	61.6 ± 8.0
	LYM# × 10 ⁹ /L	1.25 ± 1.0	2.19 ± 1.5	3.30 ± 2.0
	LYM %	36.1 ± 10.0	27.0 ± 9.0	17.8 ± 8.0
	MON# × 10 ⁹ /L	0.24 ± 0.20	0.62 ± 0.50	1.11 ± 1.10
	MON %	7.0 ± 4.0	7.7 ± 5.0	6.0 ± 6.0
	EOS# × 10 ⁹ /L	0.39 ± 0.35	0.98 ± 0.80	2.52 ± 1.30
	EOS %	11.3 ± 5.0	12.1 ± 6.0	13.6 ± 7.0
	BAS# × 10 ⁹ /L	0.10 ± 0.10	0.10 ± 0.10	0.10 ± 0.10
	BAS %	1.0 ± 1.0	1.0 ± 1.0	1.0 ± 1.0
	RBC × 10 ¹² /L	2.36 ± 0.30	4.58 ± 0.40	5.19 ± 0.50
	HGB g/L	58 ± 6	130 ± 8	163 ± 8
	HCT %	16.5 ± 3.0	37.2 ± 4.0	46.7 ± 5.0
	MCV fL	69.9 ± 5.0	81.3 ± 5.0	89.9 ± 6.0
	MCH pg	24.5 ± 2.5	28.3 ± 2.5	31.4 ± 2.5
	MCHC g/L	350 ± 30	348 ± 30	350 ± 30
	RDW-CV %	15.9 ± 3.0	14.7 ± 3.0	13.9 ± 3.0
	RDW-SD fL	39 ± 10	41 ± 10	42 ± 12
	PLT × 10 ⁹ /L	69 ± 30	273 ± 45	498 ± 70
MPV fL	10.8 ± 3.0	10.0 ± 3.0	10.0 ± 3.0	
PDW fL	10.9 ± 4.0	9.8 ± 4.0	9.7 ± 4.0	
PCT %	0.07 ± 0.07	0.27 ± 0.20	0.50 ± 0.20	
PLCR %	37.0 ± 10.0	30.5 ± 10.0	29.6 ± 10.0	
PLCC × 10 ⁹ /L	25 ± 25	83 ± 25	147 ± 40	

【NOTE】

- 1) The controls should be stored in the refrigerator (2~8°C). After opening, it will keep stable for 14 days when it is stored airtight at 2~8°C.
- 2) Please equilibrate the controls to room temperature (15~30°C) before using it.
- 3) Controls must be well mixed before use. Please mix gently to avoid cells rupturing or generating bubbles.
- 4) After using, put the controls back into the refrigerator (2~8°C) to prevent contamination and evaporation.

P/N: 01.54.459435

MPN: 01.54.459435010