


# CBC-CAL PLUS

## HEMATOLOGY CALIBRATOR

### CALIBRATOR

Assay values and expected ranges

**LOT** PLUS0225  
 2025/03/05

Method	Parameter	Calibration Values	Acceptable Range
hemaline 500 V4.01.160421 (WIC- FPGA)	WBC 10 <sup>9</sup> /L	10.33	10.13 - 10.53
	WOC 10 <sup>9</sup> /L	9.10	8.90 - 9.30
	RBC 10 <sup>12</sup> /L	4.09	4.04 - 4.14
	HGB g/dL	13.2	13.0 - 13.4
	MCV fL	86.4	85.4 - 87.4
	PLT 10 <sup>9</sup> /L	242	232 - 252
hemaline 300 application software V4.64.170701 And above (Including V4.64)	WBC 10 <sup>9</sup> /L	9.9	9.7 - 10.1
	RBC 10 <sup>12</sup> /L	4.63	4.58 - 4.68
	HGB g/dL	13.3	13.1 - 13.5
	MCV fL	87.3	86.3 - 88.3
	PLT 10 <sup>9</sup> /L	240	230 - 250

Before using, refer to the instruction sheet for mixing directions.

**MDX MEDTEK DIAGNOSTICS (NETHERLANDS) B.V.**  
 Jollemanhof 148, 1019GW Amsterdam, The Netherlands